

COMMON WORKFLOW DESCRIPTION AN EXPERIMENTAL VIEWPOINT

S. Brinckmann, F. Chen, T. Dahmen, T. Hickel, L. Huber, S. Menon,
J. Murugan, J. Neugebauer, R. Schwaiger, H. Tsybenko, O. Waseda

IMD-1: Microstructure and Properties of Materials



ALLOW FOR COMPARISON AND REUSE OF WORKFLOWS

Experiments 1

1. Prepare Al sample
2. Measure sample
3. Tensile test
4. Analyze curves
5. Calculate Modulus

ALLOW FOR COMPARISON AND REUSE OF WORKFLOWS

Experiments 1

1. Prepare Al sample
2. Measure sample
3. Tensile test
4. Analyze curves
5. Calculate Modulus

Experiments 2

1. Prepare polymer
2. Measure sample
3. Tensile test
4. Analyze curves
5. Calculate Modulus

ALLOW FOR COMPARISON AND REUSE OF WORKFLOWS

Experiments 1

1. Prepare Al sample
2. Measure sample
3. Tensile test
4. Analyze curves
5. Calculate Modulus

Experiments 2

1. Prepare polymer
2. Measure sample
3. Tensile test
4. Analyze curves
5. Calculate Modulus

Simulations

1. Build model
2. Apply forces
3. Calculate
4. Predict curves
5. Calculate Modulus

Use same workflow steps were possible

MANY WORKFLOW LANGUAGES EXIST

Dedicated workflow languages

- Common workflow language
- Snakemake
- Nextflow
- Nexus
- Kadi4Mat – json
- pyiron-workflow

...

Programming languages:

- C
- python
- bash
- ...

**TASK: PICK ONE AND USE FOR
SIMULATIONS AND EXPERIMENTS**

FORMAL WORKFLOW DESCRIPTION EXISTS

BUSINESS PROCESS MODELING

- Visual representation
- Process analysis
- Communication and documentation
- Simulation and automation

FORMAL WORKFLOW DESCRIPTION EXISTS

BUSINESS PROCESS MODELING

- Visual representation
- Process analysis
- Communication and documentation
- Simulation and automation



Two layers of information:

- Top level: short description
- Bottom level: details

GOAL: WORKFLOW STANDARDIZATION

Experiments have two layers:

- Workflow step incl. description as top layer
- Instrument has Standard Operating Procedure

Simulations have two layers:

- Function calls as top layer
- Each function is a graph node

WORKFLOW EXAMPLE

```
from common_workflow_description import Storage, Sample, step
from analysis_steps import plot_curves, calc_E
```

Header

```
wf = Workflow('Sandia Fracture Challenge 3')
proceduresLibrary = urlparse('https://...')
storage=Storage(proceduresLibrary)
```

Define workflow

```
sample = Sample('AM_NA_05')
wf.step1 = step(storage, sample, 'metallography', {})
wf.step2 = step(storage, sample, 'tensile test', {})
```

Experiments

```
...
file_name = [v for k,v in list(wf....items())...][0][1]
wf.step6 = plot_curves(file_name, 'Strain (Gauge0)', 'Engr. Stress')
wf.step7 = calc_E(file_name, 'Strain (Gauge0)', 'Engr. Stress')
```

Analysis

```
print('Output: ', '\n '.join([f"{k}: {v}" for k,v in list(wf...)]))
wf.draw().render(filename="io_demo", format="png", cleanup=True,
```

Footer

same as pyiron-workflow style

Member of the Helmholtz Association

UNIFIED WAY TO WRITE SOP: MARKDOWN+

```
# Zeiss SEM  
- Vent  
- Open chamber  
- Mount sample and close chamber  
- Use it at |voltage|20| kV  
- Use a aperture of |aperture|3|
```

Author: Steffen Brinckmann, IMD-1, FZJ

```
<!--- version:1.0 |filename|| -->
```

EXAMPLE: INSTRUCTIONS

Instructions can have:

- no forms
- only ask for a file

Common workflow description

Execute the following action with sample: AM

Metallography

```
## Grinding
• rounding all corners
• semi auto grinding: 500 -> 800 -> 1000 -> 2400 (30s) -> 4
  - significant pressure, two directions, last paper 3-4 dir
  - starting at #1000: rpm = 100
  - cleaning after each step (sample: with much water >> eth
  - Check the scratch which was made at the previous grindin

## Auto polishing
• mount sample on sample holder for auto polisher
• cleaning all clothes & polishing sample holder with water
• 3 um diamond suspension (1 spray every 1 min)
  - 30N (force reduction) / 10min 2 step (total time: 20 min

• 1 um diamond suspension (1 spray every 1 min)
  - 1st step: 15N (force reduction) / 10min / 150 rpm
  - 2nd step: 20N (force reduction) / 5min / 150 rpm
  - blue lubricant : stage 8

• cleaning after each step
  - sample: rinsing with ethanol, ultra sonic cleaning, rins
  - holder: cleaning with water

• check the surface in OM. If the surface is not good enough

## Final polishing
```

done

Common workflow description

Execute the following action with sample: AM

Confocal Image Acquisition with the Olympus

- Turn on the PC
- Turn on the LEXT OLS4000 microscope
- Start the LEXT OLS4000 software
- Select the "Imaging" mode
- Select the objective lens with the smallest setting of 5
- Press the button "Move to load position" in the "Map scre
- The objective lens will be moved upwards to avoid contact
- Place the sample on the center of the stage
- In the "Map screen" window, press the button "Move to ori
- The sample will be moved underneath the microscope under
- Move the focus knob counterclockwise to adjust the hight
- Lock the focus knob to avoid further movement
- Set the lower height limit by clicking the "Advanced sett
- Ensure that the "Imaging" mode and the "Laser Microscope"
- Search for the area of interest on the sample surface by
- Increase the magnification step-wise by clicking on the o
- Select the desired objective lens magnification from the
- After changing the objective lens, adjust height to reach
- Click on the "Laser" button. The mode will change from th
- Only the in-focus sections will be visible in the "Laser
- Click the "Focus" tab to access the focus controls
- Click the arrow buttons to raise or lower the objective l
- Set the origin height by clicking the "Advanced settings"
- Click the "Brightness" tab to change the brightness setti
- Use the "Brightness" slider or click the "Auto" button to
- Adjust the brightness settings to be as high as possible

filename

select file

done

INSTRUCTION OPTIONS

Instructions can include parameters

- geometry

- real process parameter

prescribed flow rate: 16

real flow rate: 10

- experimental metadata that is not in output file

- Setup controll box at instrument
 - General information
 - buttons F1..F3 correspond to three symbols above
 - PC-mode: use for control by PC
 - turn knob to get there and then press ESC to get move th
- Start Software and hope that LAN connection works
 - otherwise retry
- Button Bike "Travel order"
 - distance-controlled, relative: 100mm with 100mm/min
 - direction clear
- Measure the sample:
 - |width||
 - |thickness||
- Mount sample: top first and then move into bottom
- Button Setup experiment: Hammer button
 - "Test procedure, speeds", "End of test" most important
 - Switch from F0->F1 confusing
 - choose F0 as max force
 - Example speeds:
 - metals=0.5mm/min
 - polymers=5mm/min
 - Adopt path: "Store data to"
- all 3 tare-buttons should be pressed
- Button Traffic-Light: start experiments
- Button Diagram (one line)
 - open & close: rescale graph

width

1.016

thickness

1.027

filename

select file


done

EXAMPLE: ANALYSIS STEP

All possibilities exist for:

- graphical plotting (tensile curves)
- calculating material properties
- simulation steps

Young's modulus in GPa:
average=159.82
std-dev=9.05



OUTPUT AND LOGFILE

Started minimalistic workflow engine

Output:

```
step1: [Name: AM05, '', {}]  
step2: [Name: AM05, '', {'magnification': '5'}]  
step3: [Name: AM05, '....xlsx', {'width': '1.016'}]  
step4: [Name: AM05, '', {'magnification': '5'}]  
step5a: [Name: AM05, '', {'voltage': '30', 'aperture': '3'}]  
step5b: [Name: AM05, '', {'voltage': '30', 'aperture': '2'}]
```

```
08-22 10:52:39|INFO:Start workflow  
08-22 10:56:31|INFO:Start step  
    sample:{AM_NA_05}  
    procedure-name:{tensile test}  
    sha256:{ae9c351e9fcdbd41fc39d0e3831d9dc8d3388bd7e714869aa050b870dc21ac341}  
    Parameters: {  
08-22 10:57:40|INFO:Save step  
    file-name:...LongitudinalTensileOverall.xlsx  
    metadata: {"width": "1.016"}  
08-22 10:57:40|INFO:End step
```

TWO ENGINES EXIST TO EXECUTE WORKFLOWS

On desktop computer, technical staff:

- | | |
|---|-------------------------------|
| 1) writes SOPs in Word, ... | saved in Markdown+ |
| 2) define the planned workflow with GUI | saved as pyiron-workflow file |

Workflow (incl. SOP) are executed in Android App or laboratory computer

Output and log-file are saved back into ELN

FUTURE GOAL

Pyiron-workflow:

- Many features: graphics, super-computer...
- Allow for adaptive workflows

**Same output and
logging files**

Minimalistic engine:

- Can only run common-workflow-description
- No installation required
- people use reference to implement into own ELN, ...

A COMMON WORKFLOW DESCRIPTION EXISTS FOR EXPERIMENTAL AND NUMERICAL MATERIALS SCIENCE

- Most basic version exists
 - creates a receipt of all experimental steps
 - unified way to write procedures
- You can use pyiron-workflow to adopt it
- **GUI version with drag-drop will exist**

Looking for scientists feedback: What is missing?

NON UNIQUE DEFINITION: WHAT IS “A STEP” (same in simulations)

Smallest step that make scientific sense to execute individually
(might depend on others)

- SEM
- EBSD (requirement of SEM before)
- EDX (requirement of SEM before)

Issue: requirements can only be resolved if there is step types (SEM procedures, ...) → much more complicated (second version)

KNOWN SHORTCOMING

Validation of experimental input:

“What happens if I enter a non-nonsensical value?”

Issue: How to integrate that into procedure that it is easy to write?

- allowed values
- only integers can be entered
- regex needed?

KNOWN SHORTCOMING

Aborting a workflow: if sample breaks

NOT A SHORTCOMING

Allow for comment / remark in each step?

- one can enter unforeseeable events

Very easy fix: just add “|comment|” at end of file

NOT A SHORTCOMING

Do we need to have a common definition of the file? Likely yes.

But first we need a version that works, which can be tested.

Only then we find short-comings...

Development by prototyping (Wikipedia → Software prototyping)